


Configuration FACS ARIA III

Laser Bleu (488nm)


- FITC
- GFP
- Alexa Fluor 488
- CFSE
- Zs Green



- PerCP-Cy5.5
- PerCP eFluor710
- BB700


Configuration FACS ARIA III

Laser Yellow-Green (561nm)




Configuration FACS ARIA III

Laser Rouge (633nm)



Configuration FACS ARIA III

Laser Violet (405nm)



Laser	Violet - 405 nm			Blue - 488 nm			Yellow Green - 561 nm				Red - 633 nm				
	C	B	A	C	B	A	D	C	B	A	C	B	A		
BP	450/40	530/30	605/12	488/10	530/30	695/40	582/15	610/20	670/14	780/60	660/20	730/45	780/60		
Nom dans DIVA	Pacific Blue	Am Cyan	Qdot 605	SSC	FITC	PerCP-Cy5-5	PE	PE-TexasRed	PE-CY5-5	PE-Cy7	APC	Alexa Fluor 700	Alexa Fluor 770		
Autres Fluorochromes possibles	BV421	Krome Orange	BV605		GFP	PerCPeFluor 710	DsRed	IP	mNeptune	PE-Vio770	Alexa 647	Alexa Fluor 700	Alexa 750		
	eFluor450	Pacific Orange	eFluor 605		Alexa 488	BB700	Tomato	mCherry			DiD		eFluor 780		
	Vybrant Violet				CFSE								Cy5.5	APC-Vio770	
	Live Dead Violet				Zs green										APC-Cy7