

Configuration Attune: 4 lasers (405-488-561-638nm)

Lasers	Channel	Filter (BP)	Labels
Bleu (488nm)		488/10	SSC
	BL1	530/30	FITC, Alexa Fluor 488, GFP, CFSE, Sytox Green
	BL2	590/40	PE-TexasRed, IP
	BL3	695/40	PerCP, PerCP-Cy5.5, PE-AF700
Yellow (561nm)	YL1	585/16	RPE, Alexa Fluor 568, Sytox Orange
	YL2	620/15	PE-TexasRed, mCherry, Sytox AADvanced
	YL3	695/40	PE-Cy5.5, PE-AF700, PerCP-Cy5.5
	YL4	780/60	PE-Cy7, Qdot 800
Red (638nm)	RL1	670/14	APC, Alexa Fluor647, Sytox Red
	RL2	720/30	Alexa Fluor 700, Alexa Fluor 680
	RL3	780/60	APC-AF750, APC-Cy7, APC-H7
Violet (405nm)	VL1	440/50	BV421, Pacific Blue, Sytox Blue, Cell Trace Violet
	VL2	512/25	BV510, Pacific Green, CFP
	VL3	603/48	BV605, Pacific Orange, Qdot 605
	VL4	710/50	BV711, Qdot 700